

Matronix[®]

Auto Lensmeter

SLK 5600



LM 20 Touch



SLK-5600 & LM - 20

Auto Lensmeter Operation Interface

The SLK 5600 and LM 20 Touch has an elegant and contemporary style that is enhanced with a user-friendly interface that makes using it simple for both beginners and experts. With just a few clicks on its high-resolution touchscreen display, taking precise measurements is easy, cutting down on error margin and guaranteeing reliable readings every time.



SLK-5600



LM - 20



Specification SLK-5600

Brand	Matronix
Weight	3.6kg
Prism	0.00-15.00 (Step 0.01 Delta)
Display	5.6" Color LCD Screen
Printer High Speed	Thermal Printer
Power Supply	AC(110~220V) +10% 50/60Hz
Dimensions	194mm x 210mm x 408mm
Usage Application	Hospital
Color	White
Add	-10.00~+10.00D (Add, Ad 2) Step: 0.01D,0.06D,0.12D,0.25D)
Axis	-0 Degree ~180 Degree
Dot method	Ink pas
PD	Single eyes PD 20mm -39.5mm Double eye
Spherical Lens	25.00D~+25.00D (Reading increment. Step:0.01D,0.06D,0.12D,0.25)
UV Penetrability	1%~5% (Step: 1%, 5%)

Specification LM- 20

Brand	Matronix
Weight	3.6kg
Prism	0.00-15.00 (Step 0.01 Delta)
Monitor	4.7" multicolor monitor screen
Printer	High Speed, Thermal Printer
Power Supply	AC (110~220V) $\pm$ 10% 50/60Hz
Dimensions	194mm x 210mm x 408mm
Application lens diamete Add	20mm100mm-10.00+10.00 (Add, Ad 2)
Axis	Step: 0.01D, 0.06D, 0.12D, 0.25D) 0 Degree ~ 180
Dot method	Degree step 1 Degree
PD	Ink pass Single eyes PD 20mm~39.5mm
Spherical Lens	Double eye PD: 40mm85mm-25.000+25.00D (Reading increment: Step:0.01D, 0.06D, 0.12D, 0.250)
UV Penetrability	1%~5% (Step: 1%, 5%)
Cylinder lens axis	0°~180° (Reading increment step:1°

Delhi Branch (Head Office)

Kolkata Branch

Mumbai Branch

 www.matronix.in
 office@matronix.in
 01141436108
 +91 9654910394, +91 8700984670
 A-14, Okhla Industrial Area Phase-1,
Near Crowne Plaza Hotel, New Delhi-110020

 kolkata@matronix.in
 8700984032
 Chinar Park, New Town,
Kolkata-700 136

 mumbai@matronix.in
 8700984136
 1B-100, 1st Floor, Phoenix
Paragon Plaza, Kurla (W),
Mumbai 400 070