

Matronix®

SLIT LAMP

2 Step




Ziess




Haag

SLIT LAMP 2 Step

Haag



Ziess




Ziess & Haag Model

Features

- Slit Image Rotation
 - Tilting Illumination up to
 - Filter Disc
 - Halogen Lamp
 - Slit Image Width
 - Slit Image Height
 - Power Supply
- 0° to 180°
 - 20°
 - Cobalt Filter, Green Filter, Neutral
 - 12 Volt 30 Watts (Optional LED)
 - 0 to 14 mm continuously variable
 - 0 to 14 mm continuously variable
 - ~110-230V, 50-60Hz

Standard Accessories



Specification

• Magnification	2-step (10x & 16x)
• Eyepiece	12.5x wide field
• Light Source	Halogen or LED (both options available)
• Filters	Cobalt Blue, Red-Free, Heat Absorbing
• Slit Width/Length	0-14 mm adjustable
• Movement	Smooth joystick control
• Power	100-240V, 50/60Hz
• Different aperture size adjustment	Available

- | | |
|-------------------------------|--------------|
| • Eye Pieces | • 10x, 10x |
| • Objective | • 121.0 |
| • working Distance | • 100 mm |
| • Objective Lens Focal Length | • F 65 mm |
| • Real Field of View | • 20.12.5 mm |
| • Diopter Adjustment Range | • -6D to +6D |
| • Inter Pupillary Distance | • 49-90 mm |

Delhi Branch (Head Office)

 www.matronix.in
 office@matronix.in
 01141436108
 +91 9654910394, +91 8700984670
 A-14, Okhla Industrial Area Phase-1,
 Near Crowne Plaza Hotel, New Delhi-110020

Kolkata Branch

 kolkata@matronix.in
 8700984032
 Chinara Park, New Town,
 Kolkata-700 136

Mumbai Branch

 mumbai@matronix.in
 8700984136
 1B-100, 1st Floor, Phoenix
 Paragon Plaza, Kurla (W),
 Mumbai 400 070