



HAAG



ZIESS


SLIT LAMP 3 STEP




SLIT LAMP 3 STEP



HAAG



ZIESS





Ziess & Haag Model



Key Features

3-Step Magnification

10x, 16x, 25x, 40x (with 10x eyepiece)

Dual Illumination

LED & Halogen light options

Integrated Filters

Blue, Halogen, Green (Red-Free), and Diffuser

Adjustable Slit Width

0 – 12 mm (continuous)

Adjustable Slit Length

1 – 12 mm (continuous)

Slit Rotation

Full 0°–180° range

Slit Inclination Angles

0°, 5°, 10°, 15°, 20°

Universal Power Supply

100–240 V AC, 50/60 Hz

Compact Design

Only 12 kg including headrest



ACCESSORIES



1

Testing Rod

2

Fuse

3

Power Supply



Specification



• Magnification	3-step (10x & 16x)
• Eyepiece	12.5x wide field
• Light Source	Halogen or LED (both options available)
• Filters	Cobalt Blue, Red-Free, Heat Absorbing
• Slit Width/Length	0-14 mm adjustable
• Movement	Smooth joystick control
• Power	100-240V, 50/60Hz
• Different aperture size adjustment	Available

- | | |
|-------------------------------|--------------|
| • Eye Pieces | • 10x, 10x |
| • Objective | • 121.0 |
| • working Distance | • 100 mm |
| • Objective Lens Focal Length | • F 65 mm |
| • Real Field of View | • 20.12.5 mm |
| • Diopter Adjustment Range | • -6D to +6D |
| • Inter Pupillary Distance | • 49-90 mm |

Delhi Branch (Head Office)

 www.matronix.in
 office@matronix.in
 01141436108
 +91 9654910394, +91 8700984670
 A-14, Okhla Industrial Area Phase-1,

Kolkata Branch

 kolkata@matronix.in
 8700984032
 Chinar Park, New Town,
Kolkata-700 136

Mumbai Branch

 mumbai@matronix.in
 8700984136
 1B-100, 1st Floor, Phoenix
Paragon Plaza, Kurla (W),
Mumbai 400 070