



Ziess




Haag


5 STEP SLIT LAMP

SLIT LAMP 5 Step

Haag



Ziess



Ziess & Haag Model

Features

- Slit Image Rotation
 - Tilting Illumination up to
 - Filter Disc
 - Halogen Lamp
 - Slit Image Width
 - Slit Image Height
 - Power Supply
- 0° to 180°
 - 20°
 - Cobalt Filter, Green Filter, Neutral
 - 12 Volt 30 Watts (Optional LED)
 - 0 to 14 mm continuously variable
 - 0 to 14 mm continuously variable
 - ~110-230V, 50-60Hz

Standard Accessories

1

Testing Rod

Fuse

2

3

Power Supply

Specification

• Magnification	5-step (10x & 16x)
• Eyepiece	12.5x wide field
• Light Source	Halogen or LED (both options available)
• Filters	Cobalt Blue, Red-Free, Heat Absorbing
• Slit Width/Length	0-14 mm adjustable
• Movement	Smooth joystick control
• Power	100-240V, 50/60Hz
• Different aperture size adjustment	Available

- | | |
|-------------------------------|---------------------------|
| • Eye Pieces | • 10x, 10x, 16x, 25x, 40x |
| • Objective | • 121.0 |
| • working Distance | • 100 mm |
| • Objective Lens Focal Length | • F 65 mm |
| • Real Field of View | • 20.12.5 mm |
| • Diopter Adjustment Range | • -6D to +6D |
| • Inter Pupillary Distance | • 49-90 mm |

Delhi Branch (Head Office)

🌐 www.matronix.in
 ✉ office@matronix.in
 ☎ 01141436108
 📞 +91 9654910394, +91 8700984670
 📍 A-14, Okhla Industrial Area Phase-1,
 Near Crowne Plaza Hotel, New Delhi-110020

Kolkata Branch

✉ kolkata@matronix.in
 ☎ 8700984032
 📍 Chinara Park, New Town,
 Kolkata-700 136

Mumbai Branch

✉ mumbai@matronix.in
 ☎ 8700984136
 📍 1B-100, 1st Floor, Phoenix
 Paragon Plaza, Kurla (W),
 Mumbai 400 070